

Cure Ring Optical Accessory




360° of Curing Power

Cure Ring technology allows a light guide 360° of curing power. The standard ring for use with liquid light guides is available in solid or slotted versions. The hinged Cure Ring coupled to a special dual-leg fused-silica light guide easily opens and closes using a spring-loaded hinge design for rapid insertion and removal of parts. Cure Rings are ideal for bonding a number of medical devices such as catheters and tubes, cables or any other parts that require 360° curing.

Available in solid, slotted or hinged versions

Provides a 360-degree beam of light

Couples to standard OmniCure® Liquid Light Guides and High Power Fiber Light Guides



OmniCure Cure Ring Optical Accessory

Applications

This optical accessory is ideal for applications including:

- Connectorizing fibers
- Tubing or cable to connector bonding
- Catheter bonding

SPECIFICATIONS

DESCRIPTION

Fiber	Randomized fiber bundle split into ten nodes for the solid Cure Ring and nine nodes for the slotted version (Numerical Aperture = 0.22)
Throughput	Approximately 30% of the output of a single liquid filled light guide
Uniformity	+/- 10% across the nodes


PART NUMBER

PRODUCT DESCRIPTION

810-00048	Cure ring, slotted, 63.5 mm OD, 10 mm ID for 5 mm LG Typical delivery time: 10 to 12 weeks
810-00049	Cure ring, no-slot, 63.5 mm OD, 10 mm ID for 8 mm LG Typical delivery time: 10 to 12 weeks
810-00050	Cure ring, no slot, 63.5 mm OD, 15 mm ID for 5 mm LG Typical delivery time: 10 to 12 weeks

CURE RING DIMENSIONS:

All dimensions are shown in inches and millimeters. Dimensions relate to the standard OmniCure Cure Ring and are for the purpose of providing fixturing and clearance detail. Customized solutions are available upon request.



For Scandinavia and the Baltics sales and support,
please contact EFSEN UV & EB TECHNOLOGY.

Skovlytoften 33 | DK-2840 Holte
efsen@efsen.dk | phone: +45 45650260

EFSEN
UV & EB TECHNOLOGY

EXCELITAS
TECHNOLOGIES

www.excelitas.com
omnicure@excelitas.com

2260 Argentia Road
Mississauga, Ontario
L5N 6H7 CANADA
Telephone: +1 905 821-2600
Toll Free (USA and CAN): +1 800 668-8752
Fax: +1 905 821-2055

For a complete listing of our global offices, visit www.excelitas.com/locations.

©2015 Excelitas Canada Inc. OmniCure® is a registered trademark of Excelitas Canada Inc. The Excelitas logo and design are registered trademarks of Excelitas Technologies Corp. All other trademarks are the property of their respective owners, and neither Excelitas Technologies Corp., its affiliates or subsidiaries, or any of their respective products, are endorsed or sponsored by or affiliated in any way whatsoever with those organizations whose trademarks and/or logos may be mentioned herein for reference purposes. Excelitas Canada Inc. reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.